

*Z-5™ ceramic  
cad/cam dental blocks*



*Manufactured in USA*



## High Quality Z-5™ CAD/CAM Dental Blocks

Z-5™ CAD/CAM Blocks are high quality pre-formed ceramic dental blocks intended for customer CAD/CAM fabrication of zirconia crowns and/or bridges for dental restorations. Z-5™ CAD/CAM Blocks are available in an unlimited variety of shapes, sizes, shades and sintering options (pre-sintered, fully sintered, and hot isostatically pressed) according to your individual needs.

To ensure that dental professionals world-wide have the highest quality material, our US manufactured blocks are:

- FDA 510(k) cleared
- CE certified

C5 Medical Werks is a wholly owned subsidiary of CoorsTek, the largest manufacturer of technical ceramics in North America. Please contact us at [blocks@c5medicalwerks.com](mailto:blocks@c5medicalwerks.com) for immediate price and delivery information. Credit card orders are accepted.

## Available Sizes (Pre-Sintered)

C5 Part #	Forming Method	Qty in Box	Qty in Case	Shape	Dimensions (mm)
<b>Blocks</b>					
20172	Isopress	16	96	block	42 x 20 x 16
20173	Isopress	16	96	block	42 x 20 x 20
20174	Isopress	16	96	block	60 x 20 x 25
20201	Isopress	16	96	block	40 x 15 x 14
20202	Isopress	16	96	block	55 x 19 x 15.5
20203	Isopress	8	24	block	65 x 25 x 22
20211	Isopress	16	96	block	20 x 15 x 14
20212	Isopress	16	96	block	20 x 19 x 15.5
20213	Isopress	8	24	block	39 x 19 x 15.5
<b>Cylinders</b>					
20169	Uniaxial press	32	192	cylinder	16 dia. x 16
20170	Isopress	32	192	cylinder	16 dia. x 16
20171	Isopress	32	192	cylinder	20 dia. x 20
<b>Discs</b>					
20157	Isopress	2	12	disc	98.4 dia. x 10
20158	Isopress	2	12	disc	98.4 dia. x 14
20159	Isopress	2	12	disc	98.4 dia. x 16
20160	Isopress	2	12	disc	98.4 dia. x 18
20161	Isopress	2	12	disc	98.4 dia. x 20
20162	Isopress	1	6	disc	98.4 dia. x 25
20163	Isopress	2	12	disc	100 dia. x 10
20164	Isopress	2	12	disc	100 dia. x 12
20165	Isopress	2	12	disc	100 dia. x 14
20166	Isopress	2	12	disc	100 dia. x 16
20167	Isopress	2	12	disc	100 dia. x 18
20168	Isopress	2	12	disc	100 dia. x 20

Custom shapes and size are available. Please contact us at [blocks@c5medicalwerks.com](mailto:blocks@c5medicalwerks.com) for additional information.

### Advantages of C5 Medical Werks

- ISO 13485 certified quality system
- GMP compliant, FDA registered facility
- State-of-the-art ceramic manufacturing
- Manufactured in a “medical only” facility
- Credit card sales accepted

### Advantages of Z-5™ Ceramic

- Medical implant grade zirconia
- Biocompatible, state-of-the-art dental material
- High fracture durability
- Uniform shrinkage
- Manufactured in the USA
- Blocks available “off the shelf”
- White and shaded materials deliver excellent aesthetic results

### Sizes

- Standard sizes stocked for immediate delivery
- Custom blanks made to order



Manufactured in **USA**